



Lifting tub AVERO Classic



Lifting tub AVERO Classic

... the ergonomically shaped care tub

The lifting tub AVERO Classics is an ergonomically shaped care tub made out of glass-fibre-polyester. It is freely accessible from three sides and is suitable for the usage with nursing lifts. The shape of the tub body offers a safe but nevertheless comfortable ergonomic design with its integrated arm rests and shaped back rest. The safety handles inside the tub provide assistance for the entrance and exit of the patient. In the lower leg area the tub body has inserted grooves to accommodate a tub shortener.

The tub is completely closed via external paneling and is designed without any desk construction. The lifting, which is effected by two synchronized electromechanical actuators, ranges from 680 to 1.080 mm and is regulated by a cable-remote-control. A set of adjustable feet is supplied with the actuators to achieve an optimal working height.

Technical data:

Lifting tub AVERO Classic

With thermostatic mixer 1/2", incl. Hand shower
 Measures L: 1.700/2.040 mm, W: 540/800 mm,
 Tub height inside: 400 mm, Lifting height 680 – 1.080 mm,
 Power supply 230 V, 50 Hz, max. 1 kW (110 V available)
 incl. drainage and overflow fittings

Ref No.: 910206550

Accessories:

Tub shortener
 Automatic tub filling
 Sound-music-system 2RC
 Colored-light-therapy system 4
 Filling fittings 1 1/2 "
 Cleaning system
 Air-Jet-Whirl device

Ref No.: 910145600

Ref No.: 910145700

Ref No.: 910609500

Ref No.: 910609600

Ref No.: 910145400

Ref No.: 910329000

Ref No.: 910324800

Am Rübenmorgen 3
 D-35582 Wetzlar-Dutenhofen
 Phone: ++49 641/92220-0
 Fax: ++49 641/92220-20
www.beka-hospitec.de

beka
 hospitec