



Sling lift CARLO Classic



Sling lift CARLO Classic

... the multi-purpose shower and transfer lift 205 kg

Multi-purpose transfer lift CARLO Classic for the safe and comfortable transfer of residents and patients. Due to stainless steel engineering also applicable as shower lift. The extremely wide lifting range ensures the lifting from the floor and the transfer to the bed without difficulty. The four point suspension for loop slings enables a comfortable and safe transfer.

- Unrestricted usage in the shower
- Load 205 kg
- Carrier frame mechanism stainless steel and hot-dip galvanized
- Four point suspension made out of stainless steel, additionally powder-coated
- Standard electrical spreadable carrier frame
- Ergonomical-shaped handle
- Double castors for optimum moving quality, two of them lockable
- Extrem lifting range
- Option: calibratable scale

Technical data:

Carrier frame: 1.300 mm x 670 mm
 Lifting range: 1.250 mm
 Patient weight: 205 Kg

Sling lift CARLO Classic	Order No. 923001000
Adjustment and fixation of suspension by electric motor	
Sling lift CARLO Classic with scale	Order No. 923002000
Adjustment and fixation of suspension by electric motor	

Accessories:

Spare storage battery	Order No. 921070600
External battery charger	Order No. 921070700

Transfer slings, bathing slings, lavatory slings in all sizes available

Am Rübenmorgen 3
 D-35582 Wetzlar-Dutenhofen
 Fon +49 (0) 641/92220-0
 Fax +49 (0) 641/92220-20
www.beka-hospitec.de

beka
 hospitec