



Stand-up/raising aid NORA



Stand-up/raising aid NORA

... the effective relief

The transfer, for instance from the bed to the wheel chair or from the wheel chair to the bathroom, seriously strains the locomotor and skeletal system of the nursing staff and is one of their most exhausting tasks. Our new stand-up/raising aid NORA offers a considerable and effective relief of the day-to-day care.

- Easy handling
- Lifting sling with anti-slip coating
- Safety belts for the legs
- Ergonomically shaped handles
- Suitable for use in the shower
- Extra soft knee pads
- Electrically spreadable carrier frame

Technical data:

Stand-up/raising aid NORA

Max. loading capacity 200 kg

Measures L: 900 mm, W: 650 – 710 mm, H: 1.140 – 1.750 mm

Storage battery 24 V

Ref. No.: 921070000

Accessories:

Stand-up/raising sling NORA size M

Stand-up/raising sling NORA size L

Stand-up/raising sling NORA size XL

Double knee pad NORA

Spare storage battery

External battery charger 24 V

Ref. No.: 921070100

Ref. No.: 921070200

Ref. No.: 921070300

Ref. No.: 921070400

Ref. No.: 930508500

Ref. No.: 930508600

Am Rübenmorgen 3
D-35582 Wetzlar-Dutenhofen
Phone: ++49 641/92220-0
Fax: ++49 641/92220-20
www.beka-hospitec.de

beka
hospitec