



Shower trolley SINA Comfort



Shower trolley SINA Comfort

... the ideal device for transport and shower

The new shower trolley SINA Comfort, completely made of stainless steel, is the ideal device for working with high-level care residents. The sturdy safety-hinges of the side rails and the extremely resilient stretcher turn this shower trolley into an indispensable aid for the daily care and safe transport of your patients.

A special scissors-type lifting system guarantees a smooth and truly vertical electric or hydraulic height regulation. Thereby the shower trolley SINA can be adjusted to the respective working height of the nursing staff and at the same time ensures an easy transfer of the patient from the bed or the wheel chair. Eight double-castors, four of them lockable, with a special pivot point bearing provide for a smooth moving of the lift even on rough surfaces.

- **Maximum stability**
- **Extra wide height adjustment by scissors-type lifting system**
- **Completely made of powder-coated stainless steel**
- **Double-castors for smooth and safe transport**
- **Stretcher positive/negative adjustable, e.g. also applicable for Trendelenburg position**
- **Extension of the stretcher by adjustable rail**

Technical data:

Max. body weight 150 kg

Shower trolley SINA Comfort HYD-200

Hydraulic height adjustment from both sides of the trolley
Measures L: 2.050 mm, W: 760 mm, H: 670 - 1.030 mm

Ref. No.: 930507100

Shower trolley SINA Comfort HYD-160

Hydraulic height adjustment from both sides of the trolley
Measures L: 1.740 mm, W: 760 mm, H: 670-1.030 mm

Ref. No.: 930507200

Shower trolley SINA Comfort EL-200

Electric height adjustment
Including hand switch and battery charger
Measures L: 2.050 mm, W: 760 mm, H: 670 - 1.030 mm

Ref. No.: 930507300

Shower trolley SINA Comfort EL-160

Electric height adjustment
Including hand switch and battery charger
Measures L: 1.740 mm, W: 760 mm, H: 670-1.030 mm

Ref. No.: 930507400

Accessories:

SINA wedge-shaped cushion

Ref. No.: 930508320

Spare storage battery 24 V

Ref. No.: 930508500

External battery charger

Ref. No.: 930508600



Completely made of stainless steel



Sturdy safety hinges



Handle for easy positive/negative adjustment of the stretcher



Extension of the stretcher by adjustable rails



Double-castors for easy and safe transfer



Flexible water drainage hose

Am Rübenmorgen 3
D-35582 Wetzlar-Dutenhofen
Phone: ++49 641/92220-0
Fax: ++49 641/92220-20
www.beka-hospitec.de

beka
hospitec