



Tilting tub INVITA



## Tilting tub INVITA

### ... the trend-setting shower and bath solution

The ideal shower and bath tub for limited mobile patients. The special shape of this shower and bath tub ensures a non-slippery access by the patient as well as the safe usage of a patient lift. The design of the tub interior and its shaped seat with a special hygienic opening facilitates an optimal intimate care. The generously created upper body area makes nursing easy and at the same time gives the patient enough freedom of movement.

Operated by an electric motor the tub can be steplessly tilted horizontally which makes it suitable for a full-body bath at an optimal working height. The tub is placed on a carrier frame with adjustable base feet. Its extremely low water consumption and the space-saving dimensions also offer a higher economic efficiency. The trendsetting tilting tub INVITA with its seating and lying functions opens new possibilities in the field of contemporary care.

- Non-slippery independent access and exit
- Safe shower possibility with option for full-body/foot bath
- No doors, applicable for lifts
- Large hygienic opening for optimal intimate care
- Low water consumption (80 – 120 liter)
- Ergonomically shaped back part with special leg positioning
- Minor space requirement

#### Technical data:

##### Tilting tub INVITA

With thermostatic mixer for tub filling and hand shower, positioned on tub rim  
 Measures: L: 1.550 mm, W: 840 mm, H: 895/1.247 mm  
 Power supply: 230 V, 50 Hz, max. 1 kW (110 V available)

Ref. No.: 910150800

##### Tilting tub INVITA

For pre-installed fittings  
 Measures: L: 1.550 mm, W: 840 mm, H: 895/1.247 mm  
 Power supply: 230 V, 50 Hz, max. 1 kW (110 V available)

Ref. No.: 910150700

#### Accessories:

Air-Jet-Whirl device

Ref. No.: 910150830

Am Rübenmorgen 3  
 D-35582 Wetzlar-Dutenhofen  
 Phone: ++49 641/92220-0  
 Fax: ++49 641/92220-20  
[www.beka-hospitec.de](http://www.beka-hospitec.de)

**beka**  
 hospitec